

5.2.1 Using the Share Screen

Share screen features allow you to add sites and invite users to access your study.


You can access the **Share** screen from either the **My Studies** screen, **Study Designer**, or the **Settings** screen.

The **Share** button is on the bottom of the study card on the **My Studies** screen. Select **Test** or **Production** to share your study.



The screenshot shows a study card for "New Cancer Drug". At the top, the title "New Cancer Drug" is displayed. Below the title, there are two status indicators: a green circle followed by "AVAILABLE in test" and a grey circle with a slash followed by "UNPUBLISHED in production". The main body of the card contains a description: "This study tests a newly developed drug that is intended to shrink malignant tumors." Below the description, the following details are listed: "Study ID: NewCancerDrug", "Study Type: Interventional, Phase I", "Expected Enrollment: 100", and "Expected Duration: 26-Aug-2019 to 30-Aug-2020". At the bottom of the card, there is a navigation bar with three buttons: a settings gear icon, a "Design" button, and a "Share" button. The "Share" button is highlighted with a red rectangular border. Below the navigation bar is a "Go" button with a dropdown arrow.

The **Share** button is also in the header bar in **Study Designer**. Select **Test** or **Production** to share your study.



Functional approval by Kate Lambert. Signed on 2024-11-06 4:37PM

Approved for publication by Paul Bowen. Signed on 2025-02-10 1:29PM

Not valid unless obtained from the OpenClinica document management system on the day of use.