

5.2.1 Using the Share Screen

Share screen features allow you to add sites and invite users to access your study.


You can access the **Share** screen from either the **My Studies** screen, **Study Designer**, or the **Settings** screen.

The **Share** button is on the bottom of the study card on the **My Studies** screen. Select **Test** or **Production** to share your study.



The screenshot shows a study card for "New Cancer Drug". At the top, it has a title "New Cancer Drug" and two status options: "AVAILABLE in test" (selected with a green dot) and "UNPUBLISHED in production" (unselected with a grey dot). Below this, there is a description: "This study tests a newly developed drug that is intended to shrink malignant tumors." followed by study details: "Study ID: NewCancerDrug", "Study Type: Interventional, Phase I", "Expected Enrollment: 100", and "Expected Duration: 26-Aug-2019 to 30-Aug-2020". At the bottom of the card, there are three buttons: a settings gear icon, a "Design" button, and a "Share" button. The "Share" button is highlighted with a red rectangular box. Below these buttons is a "Go" button with a dropdown arrow.

The **Share** button is also in the header bar in **Study Designer**. Select **Test** or **Production** to share your study.



Approved for publication by Kerry Tamm. Signed on 2020-12-15 9:02AM

Not valid unless obtained from the OpenClinica document management system on the day of use.