



5.2.3 PDATE Datatype (Month and Year Provided)

Values collected in prior version of OpenClinica would be stored in the format specified by the users locale. The default format in the **user interface** for DATE and PDATE in OpenClinica is DD-MMM-YYYY. The default format in the user interface for DATE and PDATE in OpenClinica is DD-MMM-YYYY. Therefore, when only a month and year are provided in a PDATE field, the value in the database would be MMM-YYYY. This is the only format that will be updated automatically by a liquibase script when OpenClinica is upgraded to Amethyst. A separate script will need to be created to identify and update values in a different format.

For all partial dates stored in MMM-YYYY format, the liquibase script will convert the values YYYY-MM in order to conform to the ISO 8601 standard date format.

This page is not approved for publication.