



## 2.7.3 Remove and Restore Study Events for a Subject

You can remove a Study Event if you want to prevent the information in it from being changed, but still want to access the information and potentially restore it for use in the future. You cannot remove a Study Event when its Event status is "not scheduled." To remove a Study Event for a Subject:

1. In the Subject matrix, click the View icon in the Actions column. The View Subject page opens.
2. Click the Remove icon  for the Study Event. The Remove Event from Study page opens.
3. Click Remove Event from Study.
4. Press OK in the confirmation dialog box that displays. The View Subject page opens, reflecting the Study Event's removal for that Subject.

After you remove a Study Event for a Subject:

- The original status of the Study Event (for example, Signed or Completed) will still be reflected in the user interface. For example, the View Subject page and Printable HTML Casebook will list the Study Event as Signed, Completed, etc., but the audit log will list the Removal event and the actions available for the Study Event will change after it has been removed.
- The status for all CRFs in the Study Event that have been started will be "invalid" .
- Data in CRFs for the Study Event cannot be changed.
- Discrepancy Notes can be created for Items in CRFs for the Study Event.

You can restore a Study Event that has been removed for a Subject:

1. In the Subject matrix, click the View icon in the Actions column. The View Subject page opens.
2. Click the Restore icon  for the Study Event you want to restore. The Restore Event from Study page opens.
3. Click Restore Event to Study.
4. Press OK in the confirmation dialog box. The View Subject page opens, reflecting that the Study Event was restored for that Subject.

- After restoring a Study Event, the status of the CRFs for that Study Event for the Subject are no longer "invalid."

Approved for publication by Paul Bowen. Signed on 2025-07-23 1:58PM

Not valid unless obtained from the OpenClinica document management system on the day of use.