

## **5.4.1 About CRFs in OpenClinica**

A good way to familiarize yourself with a defined CRF is to look at its existing CRF .xls file and at the way OpenClinica presents it in the web interface:

- 1. <u>Click here</u> to download a sample CRF defined for a Physical Examination, SAMPLE\_PHYSICALEXAM\_ENGLISH.xls. Save the file to your computer.
- 2. Upload the defined CRF to your OpenClinica system.
  - 1. Set the current Study to the Study level.
  - 2. From the Tasks menu, in the Monitor and Manage Data module, select CRFs. The Manage Case Report Forms page opens.
  - 3. Click Create a New CRF. The Create a New Case Report Form page opens.
  - 4. In the MS Excel File to Upload field, click Browse, then navigate to the file you downloaded in step 1, SAMPLE\_PHYSICALEXAM\_ENGLISH.xls. Select the file, then click Open.

The path to the filename is shown in the MS Excel File to Upload field.

5. Click Preview CRF Version.

The Check CRF Version Data page opens, showing the CRF in the web interface format. 6. Click Continue.

- The Create a New CRF Version Data Committed Successfully page opens.
- 3. View the defined CRF in OpenClinica:
  - 1. From the Tasks menu, in the Monitor and Manage Data module, select CRFs. The Manage Case Report Forms page opens.
  - 2. For Sample Physical Exam English Version, click the View icon. The defined CRF displays in the web interface format.
- 4. In Excel, open the CRF file you downloaded, SAMPLE\_PHYSICALEXAM\_ENGLISH.xls.
- 5. Compare the CRF when viewed in the web interface format to the defined CRF when opened in Excel.

A Sample Defined CRF when Viewed in the Web Interface:

Sam	ple Physical Exam	English				
• C	CRF Header Info					
Click the	flag icon next to an input to enter,	view discrepancy notes.	Pleas	se note that you can only :	save the notes if CRF data entr	y has already
Ex	cit ]					
	Basic (0/9) II Body(0/35)	III Oth(0/3)	•	Select to Jump	·	
Title:	Basic Information					
Visit In	formation:					
1	Date of Physical Examination:	<b></b>	2	Time of Physical Examination:		(HH:MM)
Physica	al Exam Information:					
3	Height:	(in)	4	Weight:		(lb)
5	Temperature:	(F)	6	Pulse Rate:		(per min)
7	Respiration Rate:	(per	min)			
Blood p	ressure:					

(mm)

Note that when the CRF is previewed after upload, or viewed through 'Manage CRFs', any hidden fields in the CRF are visible.

Diastolic:

(Hg)

The Same CRF when Opened in Excel, Showing the Items Tab:

8

Systolic:

۲	00		$\bigcirc$								
🖺 🛅 🗟 🔚 🚔 🔏 🖆 🛍 🎸 🏹 • 🔍 • 💆 • 🔍 • Search in Sheet 🛛 >>											
- 1	Home Layout Ta	bles Charts Smart	Art Formulas Data Rev	view 🗸 🗸 🗸	\$.₹						
P25 ‡ 🛞 ⊘ (° fx											
	A	В	С	D							
1	ITEM_NAME	DESCRIPTION_LABEL	LEFT_ITEM_TEXT	UNITS RIGHT_ITEM	TE						
2	PEDAT	Date of Physical Exam	Physical Exam Date of Physical Examination:								
3	PETIM	Time of Physical Exam Time of Physical Examination:		HH:MM							
4	HEIGHT	Height	Height:	in	U.						
5	WEIGHT	Weight	Weight:	lb							
6	TEMPERATURE	Temperature	Temperature:	F							
7	PULSE	Pulse Rate	Pulse Rate:	per min	- 11						
8	RESPIRATION	Respiration Rate	Respiration Rate:	per min	- 11						
9	SYSTOLIC	Systolic	Systolic:	mm	-11						
10	DIASTOLIC	Diastolic	Diastolic:	Hg	-11						
11	BMI	Body Mass Index	Body Mass Index:		-11						
12	APPEARANCE	Appearance	Appearance		-11						
13	APPEARANCE_COMMENTS	Appearance Comments	Comments:(Required if Abnormal)		-11						
14 15	SKIN	Skin Skin Comments	Skin		-11						
16	SKIN_COMMENTS HEENT	H/E/E/N/T	Comments:(Required if Abnormal)		-84						
10	HEENT COMMENTS	H/E/E/N/T Comments	Comments:(Required if Abnormal)		-14						
17					1						
		Sections Groups Items	s Instructions 🗸 🗨 🔵	)	•						
	Normal View Rea	dy			11.						

This page is not approved for publication.